




Synthecon

Company Overview

Synthecon Sutures Manufacturing SA cc is a leading manufacturer of quality surgical sutures. Our sutures are recommended by top surgeons and clinicians in Southern, Central, East and West Africa. Synthecon Sutures can provide you with off the shelf or customer specific products. Through the use of only the best needles and thread in manufacturing our sutures, we meet customer requirements and guarantee the safety of patients.

www.synthecon.co.za


SYNTHECON SYNCRYL PGLA



Name: Synthecon Syncryl PGLA
Suture Type: Absorbable
Suture Lengths: 30 – 180 cm
Suture Colours: Violet, Undyed
USP Sizes: 2 – 10/0
Needle Specs: 1/2, 1/4, 3/8, 5/8, J-Needle, Straight
Fast Moving Codes: V01-25, V01-18, V01-32, V20-10, V02-32, V01-23




SYNTHECON RAPIDE



Name: Synthecon Rapide
Suture Type: Absorbable
Suture Lengths: 30 – 100 cm
Suture Colours: Violet, Undyed
USP Sizes: 1 – 8/0
Needle Specs: 1/2, 3/8, 5/8, Straight
Fast Moving Codes: 043-05, 023-05, 023-14, 023-19, 033-07, 053-45




SYNTHECON PLAIN CATGUT



Name: Synthecon Plain Catgut
Suture Type: Absorbable
Suture Lengths: 45 – 150 cm
Suture Colours: Clear
USP Sizes: 1 – 6/0
Needle Specs: 1/2, 3/8, Straight
Fast Moving Codes: 100-32, 101-32, 100-24, 130-05, 140-05, 101-14




SYNTHECON CHROMIC CATGUT



Name: Synthecon Chromic Catgut
Suture Type: Absorbable
Suture Lengths: 80 – 180 cm
Suture Colours: Brown
USP Sizes: 2 – 7/0
Needle Specs: 1/2, 1/4, 3/8, 5/8, Straight
Fast Moving Codes: 201-18, 201-42, 202-42, 201-23, 201-32, 220-41




SYNTHECON SYNTHABS PGA




Name: Synthecon Synthabs PGA
Suture Type: Absorbable
Suture Lengths: 10 – 240 cm
Suture Colours: Violet, Undyed
USP Sizes: 2 – 10/0
Needle Specs: 1/2, 1/4, 3/8, 5/8, J-Needle, Straight
Fast Moving Codes: 301-32, 301-23, 320-12, 320-10, 302-32, 330-17




SYNTHECON PDO



Name: Synthecon PDO
Suture Type: Absorbable
Suture Lengths: 45 – 240 cm
Suture Colours: Violet
USP Sizes: 2 – 6/0
Needle Specs: 1/2, 3/8, J-Needle, Straight
Fast Moving Codes: 401-60, 430-10, 420-10, 401-32, 401-23, 430-06




SYNTHECON SILK




Name: Synthecon Silk
Suture Type: Non-Absorbable
Suture Lengths: 30 – 180 cm
Suture Colours: Black, White
USP Sizes: 4 – 8/0
Needle Specs: 1/2, 1/4, 3/8, Straight
Fast Moving Codes: 520-17, 520-73, 520-10, 520-19, 520-12, 560-11



SYNTHECON POLYESTER




Name: Synthecon Polyester
Suture Type: Non-Absorbable
Suture Lengths: 30 – 150 cm
Suture Colours: Green, White
USP Sizes: 5 – 6/0
Needle Specs: 3/8, 1/4, 1/2, Straight, J-Needle
Fast Moving Codes: 622-10, 632-72, 632-10, 642-02, 642-04, 622-02




SYNTHECON POLYPROPYLENE



Name: Synthecon Polypropylene
Suture Type: Non-Absorbable
Suture Lengths: 15 – 100 cm
Suture Colours: Blue
USP Sizes: 2 – 10/0
Needle Specs: 3/8, 1/4, 1/2, Straight, J-Needle
Fast Moving Codes: 742-10, 752-02, 720-73, 720-10, 762-82, 742-04



SYNTHECON NYLON



Name: Synthecon Nylon
Suture Type: Non-Absorbable
Suture Lengths: 10 – 240 cm
Suture Colours: Blue, Black, Clear
USP Sizes: 2 – 10/0
Needle Specs: 1/2, 1/4, 3/8, 5/8, Bi-Curve, Straight
Fast Moving Codes: 820-73, 820-19, 840-05, 801-23, 830-71, 800-23



SYNTHECON MOCRYL



Name: Synthecon MoCryl
Suture Type: Absorbable
Suture Lengths: 45 – 90 cm
Suture Colours: Violet, Undyed
USP Sizes: 2 – 6/0
Needle Specs: 1/2, 3/8, 5/8, J-Needle, Straight
Fast Moving Codes: 933-57, 943-05, 901-18, 933-19, 901-32, 930-10



Contact Details

Head Office:

1 Power Street,

P.O. Box 26, Isando, 1600

Johannesburg, South Africa

Tel: 011 974 3220 | Fax: 011 974 5838

Email: sales@synthecon.co.za

Website: www.synthecon.co.za

SAMED 
 advancing patient care through MedTech

MDMSA

SAHPRA
 South African
 Health Products
 Regulatory Authority

SABS
 APPROVED


 MANAGEMENT SYSTEMS SOLUTIONS